




# UNIVERSAL SCIENTIFIC

6210 Campbell Drive. Madison, OH 44057-2003  
PHONE: 1-440-428-7800 or 1-440-829-8961 FAX: 1-440-428-8650  
E-MAIL: tomheckman@universalscientificinc.com

## Model 1502D

### Built-In Laboratory Dishwasher



#### Product Description

This entry level laboratory dishwasher has 3 cycles (Heavy, Normal, and Light) and 3 options (delay start, High Temp Wash, & Heated Dry).

- Have your laboratory dishwasher work for you, when you want it to - By choosing the delay start option, your laboratory dishwasher is programmed to start when you want it to with a 2, 4 or 6 hour delay.
- Easy to use and easy to push - with easy-clean touch pads on the laboratory dishwasher console, you can choose your desired cycle/option with the touch pads.
- Save time within your day - Instead of pre-rinsing your labware, let this laboratory dishwasher work for you with its 100% filtration which filters 100% of the water for a cleaner wash environment.

Save Water - Save Electricity - Save Time

#### Product Dimensions

- Height (Adjustable) 33 1/2" – 35"
- Width 24"
- Depth (Including Door) 25"
- Depth with Door Open 90° 49 1/4"

#### Features

- One Button Start Control Design
- Door Latch
- Tall Tub Interior Design
- White Interior Color
- 5 Wash Levels/ 1 Wash Speeds
- Plastic Filter
- Removable Filter Trap
- dB Level 61
- Static Drying System
- 6 foot three prong power cord included

#### Cycles

- 3 Cycles - Heavy, Normal, Light
- Cycle Indicator Light
- Control lock

#### Additional Settings

- Heat/No Heat Dry
- Hi-Temp Wash
- 2-4-6 Hour Delay Start

#### Rack System

- Standard Upper and Lower Rack PVC Rack Design/Rack Coating
- Small Parts Basket
- Small Items Cover

#### Specifications

- 4 Leveling Legs
- Integral Air Gap on Supply
- Water Usage (Gallons) 6.0
- Voltage Rating 120V/60Hz/15A
- Amps @ 120 Volts 10
- Shipping Weight (Approx.) 85 Lbs.

#### Available Colors



White



Black



Stainless