



# UNIVERSAL SCIENTIFIC, INC.

## 1502D SPECIFICATIONS



### STYLE:

- Undercounter Laboratory Dishwasher

### 1502D FEATURES:

- One Button Start touch controls. Displays Cycle Indicator Light
- Low rinse-Aid Indicator Light
- Standard PVC Coated Wire Rack
- Upper Rack and Lower Rack-Small Parts Basket, Small Items Cover
- Water Temperature 68.33° C (155° F) Standard
- Thermal and sound insulation
- Operation at 61 decibels
- Manual-fill detergent dispenser for non-foaming powder or liquid detergent
- Manual-fill rinse aid solution dispenser
- Static, Drying System
- Three factory-set cycle programs: HEAVY, NORMAL, LIGHT, and RINSE ONLY

- High Temp Wash Setting
- Heated Dry
- Delayed Start: 2, 4, or 6 hour delay
- Door Latch
- Button Lock
- Interior Color: White
- Interior Design: Tall tub
- Adjustable leveling feet
- Small parts basket
- 6-foot three prong power cord
- One year warranty on parts and labor

### MODEL REQUIREMENTS:

- Water consumption per fill 1.2 gallons
- Water flow rate 3 gallons (12 liters) per minute
- Water pressure at the washer 20-120 PSI

### MODEL CONFORMS TO:

- UL
- CAN/CSA

### ELECTRICAL REQUIREMENTS:

- 120 Volts, 60 Hz, 15 Amps

### DIMENSIONS/WEIGHT:

- Height (adjustable): 33 1/2" – 35"
- Width: 24"
- Depth: 24"
- Depth (Including Door): 25"
- Depth (Door Open 90°): 49 1/4"
- Shipping Weight: 85 lbs.