



# UNIVERSAL SCIENTIFIC

6210 Campbell Drive. Madison, OH 44057-2003  
PHONE: 1-440-428-7800 or 1-440-829-8961 FAX: 1-440-428-8650  
E-MAIL: tomheckman@universalscientificinc.com

## Model 5434L

### Undercounter Laboratory Glassware Washer



#### Product Description

This is our premier undercounter glassware washer. If you are looking for one of the quietest, easy-to-use glassware washers for your laboratory, the 5434L is for you! This washer runs at a quiet 48 Db(A) Operation. There are six cycles to choose from with temperatures reaching as high as 82° C/180° F.

This model comes with premium upper and lower baskets with adjustable tines and fold down shelves to accommodate even the largest loads. There is a small parts basket which can hold smaller labware that would otherwise need to be washed by hand.

The 5434L Undercounter Laboratory Glassware Washer is one of the most water and energy efficient glassware washers on the planet!

Save Water - Save Electricity - Save Time with the 5434L!

[www.universalscientificinc.com](http://www.universalscientificinc.com)

#### Product Dimensions

Height (Adjustable)	34 3/8" – 36"
Width	24"
Depth	22 7/8"
Depth with Door Open 90°	49 1/2"

#### Features

- Hidden Control Panel w/LCD Display
- Extra Tall Tub Interior Design
- 6 Temperatures Settings – Up to 82° C/180° F
- Quiet 48 dB(A) Operation
- 1 – 24 Hours Delayed Start
- End-of-Program Signal
- Stainless Steel Triple Filtration Self-Cleaning System
- Turbo Drying System
- Program Memory Recalls Last Program

#### Cycles

6 Cycles – Daily, Heavy, Normal, Quick, Eco, Rinse and Hold

#### Rack System

- Adjustable Shelves
- Premium Upper Rack – 9" Clearance
- Premium Lower Rack – 13 3/4" Clearance
- Adjustable Tines
- Small Parts Basket

#### Quality

- 18:9 Stainless Tank, Filter and Spray Arms
- Anti-Jam, Auto-Reversing Drain Pump
- Balanced Door
- Check Valve
- Galvanized Steel Base Pan
- PEX (cross-linked polyethylene) Water Line
- 6 foot long power cord included

#### Specifications

Stainless Steel Hidden Heating Element	
Adjustable leveling legs	
Annual Energy Consumption	231 kWh
Water Usage (Gallons)	.7 per fill
Voltage Rating	120V/60Hz/15A
Shipping Weight (Approx.)	195 Lbs.

#### Available Colors

