




5434L Planning & Design: Undercounter Glassware Washer

Unit Dimensions

	US	METRIC
HEIGHT:	34 3/8" – 36"	873 – 914.40 mm
WIDTH:	24"	610 mm
DEPTH:	22 7/8"	581 mm
DEPTH W/DOOR OPEN:	49 1/2"	1257 mm
SHIPPING WEIGHT:	195 lbs	88.45 kg

Technical Data


Electricity:	120 V, 60 Hz, 15 Amp
Water Pressure:	18 - 176 psi
Heating Element:	1200 watt
Max Loading:	1300 watt



Cutout Dimensions

	US	METRIC
HEIGHT:	34 3/8" – 36"	873 – 914.40 mm
WIDTH:	24"	610 mm
DEPTH:	24"	610 mm

The electrical and water supplies should enter through the area indicated by the shaded area. Preferably, they should come through the right side of the machine. The access hole must be round and smooth and no bigger than 1 1/2" (38 mm) in diameter.



Specifications are subject to change without notice. See installation instructions for additional details