



# UNIVERSAL SCIENTIFIC

6210 Campbell Drive. Madison, OH 44057-2003  
PHONE: 1-440-428-7800 or 1-440-829-8961 FAX: 1-440-428-8650  
E-MAIL: tomheckman@universalscientificinc.com

## Model 5434LF

### Freestanding Laboratory Glassware Washer



#### Product Description

How about a Freestanding Laboratory Glassware Washer that is one of the quietest in the world? This washer runs at a quiet 48 Db(A) Operation. There are six cycles to choose from with temperatures reaching as high as 82° C/180° F.

This model comes with premium upper and lower baskets with adjustable tines and fold down shelves to accommodate even the largest loads. There is a small parts basket which can hold smaller labware that would otherwise need to be washed by hand.

We give you that and more with the Universal Scientific 5434LF Freestanding Laboratory Glassware Washer!

Save Water - Save Electricity - Save Time

\*Washer is shipped with floor anchors to avoid the washer tipping forward

[www.universalscientificinc.com](http://www.universalscientificinc.com)

#### Product Dimensions

Height (Adjustable)	35 1/8" – 36 3/4"
Width	25 1/8"
Depth	28"
Depth with Door Open 90°	53 1/2"

#### Features

- Hidden Control Panel w/LCD Display
- Extra Tall Tub Interior Design
- 6 Temperatures Settings – Up to 82° C/180° F
- Quiet 48 dB(A) Operation
- 1 – 24 Hours Delayed Start
- End-of-Program Signal
- Stainless Steel Triple Filtration Self-Cleaning System
- Turbo Drying System
- Program Memory Recalls Last Program

#### Cycles

6 Cycles – Daily, Heavy, Normal, Quick, Eco, Rinse and Hold

#### Rack System

Adjustable Shelves  
 Premium Upper Rack – 9" Clearance  
 Premium Lower Rack – 13 3/4" Clearance  
 Adjustable Tines  
 Small Parts Basket

#### Quality

- 18:9 Stainless Tank, Filter and Spray Arms
- Anti-Jam, Auto-Reversing Drain Pump
- Balanced Door
- Check Valve
- Galvanized Steel Base Pan
- PEX (cross-linked polyethylene) Water Line
- 6 foot long power cord included

#### Specifications

Stainless Steel Hidden Heating Element  
 Adjustable leveling legs  
 Annual Energy Consumption 231 kWh  
 Water Usage (Gallons) .7 per fill  
 Voltage Rating 120V/60Hz/15A  
 Shipping Weight (Approx.) 195 Lbs.

#### Available Colors

