



UNIVERSAL SCIENTIFIC

6210 Campbell Drive, Madison, OH 44057-2003

PHONE: 1-440-428-7800 or 1-440-829-8961 FAX: 1-440-428-8650

E-MAIL: tomheckman@universalscientificinc.com

Model 5434LP

Mobile Glassware Washer



Product Description

How about a Mobile Laboratory Glassware that you can move around to different rooms in your laboratory when needed?

This washer runs at a quiet 48 dB(A) Operation, has six cycles to choose from with temperatures reaching up to 82° C/180° F. The extra tall tub tank, filter, spray arms, and hidden heating element are all made of stainless steel. You also have the option of delaying the washer to start up to 24 hours.

This model comes with premium upper and lower baskets with adjustable tines and fold down shelves to accommodate even the largest loads. There is a small parts basket which can hold smaller labware that would otherwise need to be washed by hand.

Save Water - Save Electricity - Save Time

www.universalscientificinc.com

Product Dimensions

Height (Adjustable)	39 3/8" – 41"
Width	25"
Depth	27 7/8"
Depth with Door Open 90°	51 1/2"

Features

- Hidden Control Panel w/LCD Display
- Extra Tall Tub Interior Design
- 6 Temperatures Settings – Up to 82° C/180° F
- Quiet 48 dB(A) Operation
- 1 – 24 Hours Delayed Start
- End-of-Program Signal
- Stainless Steel Triple Filtration Self-Cleaning System
- Turbo Drying System
- Program Memory Recalls Last Program

Cycles

6 Cycles – Daily, Heavy, Normal, Quick, Eco, Rinse and Hold

Rack System

- Adjustable Shelves
- Premium Upper Rack – 9" Clearance
- Premium Lower Rack – 13 3/4" Clearance
- Adjustable Tines
- Small Parts Basket

Quality

- 18:9 Stainless Tank, Filter and Spray Arms
- Anti-Jam, Auto-Reversing Pump
- 18 gauge type 304 stainless steel sides and back
- 3" diameter rubber wheels with swivel and brake on front- rigid on back
- water supply & drain connection kit to lab sink Hot Water x Cold Water fixture
- 3/4" Thick black plastic laminate top,
- Balanced Door
- Check Valve
- Galvanized Steel Base Pan
- PEX (cross-linked polyethylene) Water Line
- 12-foot long power cord included

Specifications

Stainless Steel Hidden Heating Element	
Adjustable leveling legs	
Annual Energy Consumption	231 kWh
Water Usage (Gallons)	.7 per fill
Voltage Rating	120V/60Hz/15A
Shipping Weight (Approx.)	195 Lbs.

Available Colors



White



Black



Stainless