



UNIVERSAL SCIENTIFIC

6210 Campbell Drive. Madison, OH 44057-2003
PHONE: 1-440-428-7800 or 1-440-829-8961 FAX: 1-440-428-8650
E-MAIL: tomheckman@universalscientificinc.com

Model 5434S

Undercounter Spindle Rack Glassware Washer



Product Description

Our 5434S Undercounter Spindle Rack Laboratory Glassware Washer super-cleans even the hardest to clean flasks and glassware.

Six cycle programs give the user full control in setting water temperature, wash and rinse times, and number of rinses with temperatures reaching up to 82° C (180° F). The result is thoroughly clean, dry labware, ready for use in sensitive research.

The washer is outfitted with a Premium upper basket with fold down shelves to wash your everyday bowls, beakers, and medium sized labware. Meanwhile, the bottom has 34 Spindles for washing narrow-neck labware. In addition, a premium lower basket is included and can be used instead of the Lower Spindle Rack.

Save Water - Save Electricity - Save Time

www.universalscientificinc.com

Product Dimensions

Height (Adjustable)	34 3/8" – 36"
Width	24"
Depth	26"
Depth with Door Open 90°	53 1/2"

Features

- Hidden Control Panel w/LCD Display
- Extra Tall Tub Interior Design
- 6 Temperatures Settings – Up to 82° C/180° F
- Quiet 48 dB(A) Operation
- 1 – 24 Hours Delayed Start
- End-of-Program Signal
- Stainless Steel Triple Filtration Self-Cleaning System
- Turbo Drying System
- Program Memory Recalls Last Program

Cycles

6 Cycles – Daily, Heavy, Normal, Quick, Eco, Rinse and Hold

Rack System

Adjustable Shelves

Premium Upper Basket – 9" Clearance

Premium Lower Spindle Rack –

34 customizable length spindles. 8" standard

Additional Lower Basket –

Option to switch out Spindle Rack to wash larger wide mouth labware. 13' 3/4" Clearance

Quality

- 18:9 Stainless Tank, Filter and Spray Arms
- Anti-Jam, Auto-Reversing Pump
- Balanced Door
- Check Valve
- Galvanized Steel Base Pan
- PEX (cross-linked polyethylene) Water Line
- 6 foot long power cord included

Specifications

Stainless Steel Hidden Heating Element

Adjustable leveling legs


Annual Energy Consumption 231 kWh

Water Usage (Gallons) .7 per fill


Voltage Rating 120V/60Hz/15A

Shipping Weight (Approx.) 195 Lbs.


Available Colors



White



Black



Stainless