



UNIVERSAL SCIENTIFIC

6210 Campbell Drive. Madison, OH 44057-2003
PHONE: 1-440-428-7800 or 1-440-829-8961 FAX: 1-440-428-8650
E-MAIL: tomheckman@universalscientificinc.com

Model 5434SP

Mobile Spindle Rack Glassware Washer



Product Description

How about a Mobile Spindle Rack Laboratory Glassware that you can move around to different rooms in your laboratory when needed?

This washer runs at a quiet 48 dB(A) Operation, has 6 cycles to choose from with temperatures reaching up to 82° C/180° F. The extra tall tub tank, filter, spray arms, and hidden heating element are all made of stainless steel. You also have the option of delaying the washer to start up to 24 hours.

The washer is outfitted with a Premium upper basket with fold down shelves to wash your everyday bowls, beakers, and medium sized labware. Meanwhile, the bottom has 34 Spindles for washing narrow-neck labware. In addition, a premium lower basket is included and can be used instead of the Lower Spindle Rack.

Save Water - Save Electricity - Save Time

www.universalscientificinc.com

Product Dimensions

Height (Adjustable)	39 3/8" – 41"
Width	24"
Depth	28 7/8"
Depth with Door Open 90°	55 1/2"

Features

- Hidden Control Panel w/LCD Display
- Extra Tall Tub Interior Design
- 6 Temperatures Settings – Up to 82° C/180° F
- Quiet 48 dB(A) Operation
- 1 – 24 Hours Delayed Start
- End-of-Program Signal
- Stainless Steel Triple Filtration Self-Cleaning System
- Turbo Drying System
- Program Memory Recalls Last Program

Cycles

6 Cycles – Daily, Heavy, Normal, Quick, Eco, Rinse and Hold

Rack System

Adjustable Shelves
 Premium Upper Basket – 9" Clearance
 Premium Lower Spindle Rack – 34 customizable length spindles. 8" standard
 Additional Lower Basket – Option to switch out Spindle Rack to wash larger wide mouth labware. 13' 3/4" Clearance

Quality

- 18-9 Stainless Tank, Filter and Spray Arms
- Anti-Jam, Auto-Reversing Pump
- Balanced Door
- Check Valve
- Galvanized Steel Base Pan
- PEX (cross-linked polyethylene) Water Line
- 6 foot long power cord included

Specifications

Stainless Steel Hidden Heating Element
 Adjustable leveling legs
 Annual Energy Consumption 231 kWh
 Water Usage (Gallons) .7 per fill
 Voltage Rating 120V/60Hz/15A
 Shipping Weight (Approx.) 195 Lbs.

Available Colors

