



UNIVERSAL SCIENTIFIC

6210 Campbell Drive. Madison, OH 44057-2003
PHONE: 1-440-428-7800 or 1-440-829-8961 FAX: 1-440-428-8650
E-MAIL: tomheckman@universalscientificinc.com

Model 5634ADA

American Disabilities Act (ADA)
Laboratory Glassware Washer



Product Description

FINALLY! A Laboratory Glassware Washer that is specifically designed to fully comply with the American Disabilities Act (A.D.A.)! This ADA glassware washer features 9 different cycles to choose from (Daily, Heavy, Normal, Delicate, Quick, Super Quick, Eco, Rinse & Hold, Rinse and Dry).

-Have your ADA laboratory glassware washer work for you, when you want it to - By choosing the delay start option, your ADA laboratory glassware is programmed to start when you want it to with the option to choose from 1 to 24 hour delay.

-This washer comes with 7 temperature settings as high as 82° C/180° F. Using just .7 gallons of water per fill. This model also features a Triple Filtration Self-Cleaning System, end-of-program signal, and Hidden Control Panel w/LCD Display!

Save Water - Save Electricity - Save Time

www.universalscientificinc.com

Product Dimensions

Height (Adjustable)	32 1/4" – 34"
Width	24"
Depth	22 7/8"
Depth with Door Open 90°	48 3/8"

Features

- Hidden Control Panel w/LCD Display
- Front LED Display
- Tall Tub Interior Design
- 7 Temperatures Settings – Up to 82° C/180° F
- Quiet 48 dB(A) Operation
- 1 – 24 Hours Delayed Start
- End-of-Program Signal
- Stainless Steel Triple Filtration Self-Cleaning System
- Turbo Drying System

Cycles

9 Cycles – Daily, Heavy, Normal, Delicate, Quick, Super Quick, Eco, Rinse and Hold, Rinse and Dry

Rack System

Adjustable Shelves
 Premium Upper Rack – 9" Clearance
 Lower Rack – 12" Clearance
 Adjustable Tines
 Small Parts Basket

Quality

- 18:9 Stainless Tank, Filter and Spray Arms
- Anti-Jam, Auto-Reversing Pump
- Balanced Door
- Check Valve
- Galvanized Steel Base Pan
- PEX (cross-linked polyethylene) Water Line
- 6 foot long power cord included

Specifications

Stainless Steel Hidden Heating Element
 4 Leveling Legs
 Annual Energy Consumption 231 kWh
 Water Usage (Gallons) .7 per fill
 Voltage Rating 120V/60Hz/15A
 Shipping Weight (Approx.) 185 Lbs.

Available Colors

