




# 5634ADA Planning & Design: ADA Undercounter Glassware Washer

## Unit Dimensions

|                    | US            | METRIC       |
|--------------------|---------------|--------------|
| HEIGHT:            | 32 1/4" – 34" | 819 - 864 mm |
| WIDTH:             | 24"           | 610 mm       |
| DEPTH:             | 22 7/8"       | 581 mm       |
| DEPTH W/HANDLE:    | 24 7/8"       | 632 mm       |
| DEPTH W/DOOR OPEN: | 48 3/8"       | 1229 mm      |
| WEIGHT:            | 185 lbs.      | 84 kg        |

## Technical Data

|                  |                      |
|------------------|----------------------|
| Electricity:     | 120 V, 60 Hz, 15 Amp |
| Water Pressure:  | 18 - 176 psi         |
| Heating Element: | 1200 watt            |
| Max Loading:     | 1300 watt            |




## Cutout Dimensions

|         | US            | METRIC       |
|---------|---------------|--------------|
| HEIGHT: | 32 1/4" – 34" | 819 - 864 mm |
| WIDTH:  | 24"           | 610 mm       |
| DEPTH:  | 22 7/8"       | 581 mm       |

Depth includes space for the high-loop drain hose

The electrical and water supplies should enter through the area indicated by the shaded area. Preferably, they should come through the right side of the machine. The access hole must be round and smooth and no bigger than 2" (51 mm) in diameter.



Specifications are subject to change without notice. See installation instructions for additional details