



UNIVERSAL SCIENTIFIC

6210 Campbell Drive. Madison, OH 44057-2003
PHONE: 1-440-428-7800 or 1-440-829-8961 FAX: 1-440-428-8650
E-MAIL: tomheckman@universalscientificinc.com

Model 5634ADAS

American Disabilities Act (ADA)
Undercounter Spindle Rack Washer



Product Description

This laboratory Spindle Rack Washer fully complies with the American Disabilities Act (A.D.A.)! This ADA glassware washer features 9 different cycles to choose from (Daily, Heavy, Normal, Delicate, Quick, Super Quick, Eco, Rinse & Hold, Rinse and Dry).

-The top basket is a Premium Upper Rack with a 9" clearance. The bottom rack features 34 customizable length lower spindle rack 6" standard. We also include an additional Exclusive Lower Rack 12" clearance with adjustable tines and adjustable shelves.

-Program your ADA Spindle Rack Washer to start when you want it to with the option to choose from 1 to 24 hour delayed start.

-This washer comes with 7 temperature settings as high as 82° C/180° F. Using just .7 gallons of water per fill. This model also features a Triple Filtration Self-Cleaning System, end-of-program signal, and Hidden Control Panel w/LCD Display!

www.universalscientificinc.com

Product Dimensions

Height (Adjustable)	32 1/4" – 34"
Width	24"
Depth	28 7/8"
Depth with Door Open 90°	54 3/8"

Features

- Hidden Control Panel w/LCD Display
- Front LED Display
- Tall Tub Interior Design
- 7 Temperatures Settings – Up to 82° C/180° F
- Quiet 48 dB(A) Operation
- 1 – 24 Hours Delayed Start
- End-of-Program Signal
- Stainless Steel Triple Filtration Self-Cleaning System
- Turbo Drying System

Cycles

9 Cycles – Daily, Heavy, Normal, Delicate, Quick, Super Quick, Eco, Rinse and Hold, Rinse and Dry

Rack System

- Adjustable Shelves
- Premium Upper Rack – 9" Clearance
- Lower Spindle Rack – 34 customizable length lower spindle rack 6" standard
- Additional Lower Rack – 12" Clearance, Adjustable Tines
- Small Parts Basket

Quality

- 18:9 Stainless Tank, Filter and Spray Arms
- Anti-Jam, Auto-Reversing Pump
- Balanced Door
- Check Valve
- Galvanized Steel Base Pan
- PEX (cross-linked polyethylene) Water Line
- 6 foot long power cord included

Specifications

Stainless Steel Hidden Heating Element	
4 Leveling Legs	
Annual Energy Consumption	231 kWh
Water Usage (Gallons)	.7 per fill
Voltage Rating	120V/60Hz/15A
Shipping Weight (Approx.)	185 Lbs.

Available Colors

